



KROMATON[®]
ROUSSELET CENTRIFUGATION **GROUP**

FCPC[®]

KROMATON

Since 1999 Kromaton has pioneered in developing Fast Centrifugal Partition Chromatography (FCPC[™]) systems for the separation and purification of complex mixtures of natural, synthetic or biological/biotechnological origin.

KROMATON IS THE WORLD MARKET LEADER IN CENTRIFUGAL PARTITION CHROMATOGRAPHY (CPC)

with more than 130 systems operating around the globe.

Kromaton's technological expertise is powered by Rousselet Robatel know-how in designing and manufacturing industrial centrifugal systems made of stainless steel.

This guarantees effectiveness and robustness through state-of-the-art design.

WWW.KROMATON.COM
WWW.ROUSSELET-ROBATEL.COM

MARKETS/SECTORS OF ACTIVITY

Natural compounds / ingredients of plant origin

(Reference standards, APIs, Nutraceuticals, Nutrients, Cosmetics etc.)

Synthetic compounds

(APIs, Specialty chemicals, Pesticides etc.)

Biotechnology-derived mixtures

(Biopharmaceuticals etc.)



WHAT IS FCPC ?

Fast Centrifugal Partition Chromatography (FCPC™) is a liquid/liquid purification technology for high-added-value compounds and molecules.

It is a preparative chromatography technique based on the principle that both mobile and stationary phases are liquid. This, a biphasic solvent system is utilized and the stationary phase is maintained in the column by centrifugal force, while mobile phase percolates through, thus presenting several advantages in comparison with conventional chromatography techniques :

- No solid packing material to dispose.
- High stationary phase ratio in the column.
- No sample loss due to irreversible adsorption.
- No denaturation of sensitive compounds.
- High throughput capacity.
- Moderate solvent consumption.
- Wide range of possible applications in terms of polarity.
- Easy and predictable scale-up.
- Versatile operation (elution/extrusion, normal/reverse-phase, pH-zonerefining, ion-exchange, dual mode, gradient elution, etc).

PRINCIPLE OF OPERATION


Centrifugal Partition Chromatography operation relies on a simple principle.

A biphasic solvent system is the core of this operation :

- One phase is selected to act as a stationary phase, either the lighter or the heavier. This is achieved by centrifugal force.
- The other phase is continuously percolating through, thus acting as the mobile phase. This is achieved by a pumping system.

The solutes of interest are introduced into the chromatographic column and are migrating along its length according to their distribution coefficient between the mobile and stationary phase.

The column (also called rotor) is specifically designed with twin-cells and ducts that guarantee optimal surface of contact between the two phases.



PRODUCTS

FCPC C - LABORATORY SCALE EQUIPMENT

FCPC C is a benchtop system, ideal for method development and preparative applications involving small sample quantities, low solvent consumption and faster separations.

Two interchangeable rotors can be mounted on this skid, 25 and 50ml.

DATA	FCPC C25	FCPC C50
Rotor/Column volume	25 ml	50 ml
Sample quantity	10 - 250 mg	10 - 500 mg
Flow rate	1 - 10 ml/min	1 - 10 ml/min
Separation time	5 - 30 min	
Pressure max.	80 bars / 1,160 psi	
Rotation speed max.	3,000 rpm	
Size/Weight power	300 x 300 x h 555 mm/40 kg 750W - 110 or 220 VAC	
Material of construction	Rotor: Stainless steel + PTFE. Housing: Stainless steel	



FCPC C is Kromaton's budget solution for feasibility studies and research

FCPC A - VERSATILE & SCALABLE EQUIPMENT

FCPC A is a benchtop system, ideal for method development as well as preparative applications with upscalable results and manual or automated control powered by an integrated LC station.

DATA	FCPC A25	FCPC A50	FCPC A200	FCPC A1000
Rotor/Column volume	25 ml	50 ml	200 ml	1,000 ml
Sample quantity	10-250mg	10-500mg	100mg-5g	1g-30g
Flow rate	1-10ml/min	1-10ml/min	10-50ml/min	10-80ml/min
Separation time	5 - 20 min		30 - 180 min	
Pressure max.	80 bars / 1,160 psi			
Rotation speed max.	3,000 rpm			2,000 rpm
Size/Weight power	630 x 437 x h 630 mm/115 kg 750W - 110 or 220 VAC			
Security	Emergency stop, Overpressure, Solvent Retention Tank, Cooling System			
Material of construction	Rotor: Stainless steel + PTFE. Housing: Stainless steel			
Options	Biocompatible version (titanium). Reduced cell number for higher throughput.			




FCPC A is Kromaton's most versatile system: four interchangeable rotors, 25, 50, 200 and 1,000ml, can be mounted on the same skid

Additionally it can take a **Centrifugal Partition Extractor rotor (CPE)** of 300ml total volume.

FCPC D - INDUSTRIAL PRODUCTION SYSTEM

FCPC D is Kromaton's **GMP-compatible Production system**, for large-scale purification. It is delivered with automated control powered by an integrated LC station. It can take two interchangeable rotors, of 5 and 10L volume and can be manufactured in an **ATEX version**.

DATA	FCPC D5000	FCPC D10000
Rotor/Column volume	5 liters	10 liters
Sample quantity	50-500	100-000
Flow rate	700 ml/min	1,500 ml/min
Separation time	2 to 4 hours	
Pressure max.	60 bars / 870 psi	
Rotation speed max.	1,400 rpm	
Dimensions	1,710mm x 1,160mm x 1,480mm (height)	
Poids	1,400 kg	1,600 kg
Puissance	5.5 kW - 400 V	
Material of construction	Stainless steel 316 L / PEEK / PTFE	
Housing Material	Stainless steel 316 L	
Options	FCPC version with reduced cell number for higher throughput. ATEX	



SOLVENT SYSTEMS

Non aqueous systems
for apolar applications.

Medium polarity
Ternary/Quaternary
System families (Arizona, ChMW, etc.)

Aqueous
for polar applications.

FAST CENTRIFUGAL PARTITION CHROMATOGRAPHY



KROMATON®
ROUSSELET CENTRIFUGATION GROUP

HEADQUARTERS
ROUSSELET CENTRIFUGATION
SA 45 Avenue Rhin & Danube
PAE de Marenton
07104 ANNONAY - FRANCE
Tel. : +33 / (0) 4 75 69 22 11
Fax: +33 / (0) 4 75 67 69 80
E-mail: info@rousselet-robotel.com

UK
Rousselet UK Ltd
Parkside House, 17 East Parade
HARROGATE
NORTH YORKSHIRE HG1 5LF
Tel : + 44 (0)1 423 530 093
Fax : + 44 (0)1 423 530 120
E-mail : devans@rousselet-robotel.com

DEUTSCHLAND
ARRGOS GmbH
Max-Eyth-Str. 1
D-71691 Freiberg a.N.
Tel : +49(0)7 7229 20
Fax : +0)173 757 6226
Email : nicolas.kiehl@arrgos.de

USA
Robotel Inc.
703 West Housatonic Street
PITTSFIELD - MA 01201
Tel : + 1 413 499 4818
Fax : + 1 413 499 5648
E-mail : sales@rousselet-robotel.us
Website : https://www.rousselet-robotel.us/